

Immunofluorescence (IF/IFC) Protocol

I. Materials and Reagents:

- Cell culture medium
- Chamber slide and cell culture plate
- 1X PBS (pH7.4)
- 1X PBST (0.05% Tween 20)
- Fixing solution: 4% Formalin/dd.H₂O
- Penetration solution: 0.1% Triton X-100/1xPBS
- Blocking buffer: 0.3% BSA/1xPBS
- Primary antibody and Fluorochrome-conjugated second antibody
- Staining solutions: DAPI and Phalloidin
- Mounting medium
- Glass slides and tweezers
- Fluorescence microscopy

II. Experimental Procedure:

1. Cell pretreatment:

- Wash cultured cells with PBS 2 times gently, 5 minutes per wash.
- Fix by adding 0.2ml of 4% Formalin to each well, fix cells at RT for 20 minutes, then remove fixing solution and add approximately 0.5ml of 1X PBS to cover the cells.
- Freeze cells at -20° (if proceeding, skip this step after fixing.)
- Thaw treated cells at room temperature approximately 30 minutes.
- Wash treated cells with 1X PBS 2 times, 5 minutes per wash.

2. Penetration: penetrate cells with 0.1% Triton X-100 at room temperature for 10 minutes.

3. Wash: wash penetrated cells with 1X PBS twice, each time for 5 minutes.

4. Block: block cells with 0.3ml of 0.3% BSA at room temperature for 30 minutes.

5. Primary antibody: remove the blocking solution, add primary antibody at recommended dilution into the well and incubate at 37°C for 60 minutes.

6. Wash: remove the primary antibody and wash cells with PBST 3 times, 5 minutes per wash.

7. **Secondary antibody:** add diluted secondary antibody in blocking buffer at a desired concentration into cells, incubated at 37°C for 50 minutes. Avoid light by covering with foil paper.
8. **Wash:** remove the secondary antibody, wash cells with PBST 3 times, 5 minutes per wash.
9. **Nuclei staining:** add 0.2ml diluted DAPI solution at 10ug/ml into the well, and incubate at room temperature for 10 minutes (skip this step and go to step 13 if so desired).
10. **Wash:** remove the DAPI solution, wash cells with PBST 3 times, 5 minutes per wash.
11. **Cytoplasmic actin staining:** add 0.3ml diluted Phalloidin solution into the well, incubate at 37° C for 1 hour or 4°C overnight.
12. **Wash:** remove the Phalloidin solution and wash cells two times with PBST and one time with sterilized water, 5 minutes per wash.
13. **Mount:**
 - Use a clean glass slide to mark slide/name correctly.
 - Drop proper amount of mounting medium on the glass slide.
 - Take the chamber slide out from the well using tweezers, keep cells side down to cover on mounting medium of the glass slide.
 - Remove excess medium from the edge of chamber slide, try to avoid bubbles.
14. **Visualization:** observe cells using Fluorescence microscopy.

