

## 14-3-3 Rabbit mAb


ABZ5066

### Product Overview


<b>Name</b>	14-3-3 Rabbit mAb
<b>Catalog #</b>	ABZ5066
<b>Clonality</b>	Monoclonal
<b>Accession(Primary)</b>	P31946
<b>Application Note (Approx.)</b>	WB1:500-1:1000 IHC1:100-1:200
<b>Precautions</b>	

### Target information(P31946)

<b>Synonyms</b>	
<b>Gene ID</b>	7529/7534
<b>Other Names</b>	
<b>Function</b>	
<b>Cellular location</b>	
<b>Note</b>	



Western blot analysis of 14-3-3 in HeLa lysates using 14-3-3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using 14-3-3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.