

ANTI-ACE2 ANTIBODY (ABP1050)

Rabbit Polyclonal Antibody


ABP1050

Product Overview

Name	ANTI-ACE2 ANTIBODY (ABP1050)
Catalog #	ABP1050
Clonality	Polyclonal
Accession(Primary)	Q9BYF1
Application Note (Approx.)	WB1:1000
Precautions	

Target information(Q9BYF1)

Synonyms**Gene ID** 59272**Other Names****Function****Cellular location****Note**



Sample: Human (HepG2) cell lysate. PVDF membrane was probed with Anti-ACE2 Antibody (Catalog# PP1050-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.