

# ANTI-AIFM1 (AIF) ANTIBODY (ABP1063)

*Rabbit Polyclonal Antibody*


ABP1063

## Product Overview

<b>Name</b>	ANTI-AIFM1 (AIF) ANTIBODY (ABP1063)
<b>Catalog #</b>	ABP1063
<b>Clonality</b>	Polyclonal
<b>Accession(Primary)</b>	O95831
<b>Application Note (Approx.)</b>	WB1:1000
<b>Precautions</b>	

## Target information(O95831)

**Synonyms****Gene ID** 9131**Other Names****Function****Cellular location****Note**



Sample: Human (293) cell lysate. PVDF membrane was probed with Anti-AIFM1/AIF Antibody (Catalog# PP1063-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.