

ANTI-ALKALINE PHOSPHATASE ANTIBODY (ABP1035)

Rabbit Polyclonal Antibody

ABP1035

Product Overview

Name	ANTI-ALKALINE PHOSPHATASE ANTIBODY (ABP1035)
Catalog #	ABP1035
Clonality	Polyclonal
Accession(Primary)	P05186
Other Accession	P09242
Application Note (Approx.)	WB1:2000
Precautions	

Target information(P05186)

Synonyms


Gene ID 249

Other Names

Function

Cellular location

Note



Samples: Human (HeLa) and Mouse (F9) cell lysates.

PVDF membrane was probed with Anti-Alkaline

Phosphatase Antibody (Catalog# PP1035-S1) at 1:2000

dilution, followed by HRP conjugated secondary antibody.

Experiment performed under reducing conditions.