

ANTI-ALKALINE PHOSPHATASE ANTIBODY (ABP1035)

Rabbit Polyclonal Antibody


ABP1035

Product Overview

| | |
|-----------------------------------|-------------------------------------------------|
| Name | ANTI-ALKALINE PHOSPHATASE ANTIBODY (ABP1035) |
| Catalog # | ABP1035 |
| Clonality | Polyclonal |
| Accession(Primary) | P05186 |
| Other Accession | P09242 |
| Application Note (Approx.) | WB1:2000 |
| Precautions | |

Target information(P05186)

Synonyms**Gene ID** 249**Other Names****Function****Cellular location****Note**



Samples: Human (HeLa) and Mouse (F9) cell lysates.
PVDF membrane was probed with Anti-Alkaline Phosphatase Antibody (Catalog# PP1035-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.