

# ANTI-ALPHA TUBULIN ANTIBODY (ABP1008)

*Rabbit Polyclonal Antibody*

ABP1016

## Product Overview

Name	ANTI-ALPHA TUBULIN ANTIBODY (ABP1008)
Catalog #	ABP1016
Clonality	Polyclonal
Accession(Primary)	P68363
Application Note (Approx.)	WB1:2000
Precautions	

## Target information(P68363)

### Synonyms


Gene ID 10376

### Other Names

### Function

### Cellular location

### Note



Samples: Human (293T/17, A431, HeLa, HepG2, Jurkat, Raji) cell lysates. PVDF membrane was probed with Anti-Alpha Tubulin Antibody (Catalog# PP1016-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.