

ANTI-ALPHA TUBULIN ANTIBODY (ABP1008)

Rabbit Polyclonal Antibody


ABP1016

Product Overview

| | |
|-----------------------------------|---------------------------------------|
| Name | ANTI-ALPHA TUBULIN ANTIBODY (ABP1008) |
| Catalog # | ABP1016 |
| Clonality | Polyclonal |
| Accession(Primary) | P68363 |
| Application Note (Approx.) | WB1:2000 |
| Precautions | |

Target information(P68363)

Synonyms**Gene ID** 10376**Other Names****Function****Cellular location****Note**



Samples: Human (293T/17, A431, HeLa, HepG2, Jurkat, Raji) cell lysates. PVDF membrane was probed with Anti-Alpha Tubulin Antibody (Catalog# PP1016-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.