

# ANTI-ATF4 ANTIBODY (ABP1012)

*Rabbit Polyclonal Antibody*

ABP1012

## Product Overview

Name	ANTI-ATF4 ANTIBODY (ABP1012)
Catalog #	ABP1012
Clonality	Polyclonal
Accession(Primary)	P18848
Application Note (Approx.)	WB1:1000
Precautions	

## Target information(P18848)

### Synonyms


Gene ID 468

### Other Names

### Function

### Cellular location

### Note



Samples: Human (Jurkat, MCF-7) cell lysates and Mouse (heart) and Rat (heart) tissue lysates. PVDF membrane was probed with Anti-ATF4 Antibody (Catalog# PP1012-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.