

ANTI-ATG12 ANTIBODY (ABP1039)

Rabbit Polyclonal Antibody


ABP1039

Product Overview

Name	ANTI-ATG12 ANTIBODY (ABP1039)
Catalog #	ABP1039
Clonality	Polyclonal
Accession(Primary)	O94817
Application Note (Approx.)	WB1:2000
Precautions	

Target information(O94817)

Synonyms**Gene ID** 9140**Other Names****Function****Cellular location****Note**



Samples: Human (A549, HCT1116, HeLa) cell lysates.
PVDF membrane was probed with Anti-ATG12 Antibody (Catalog# PP1039-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.