

ANTI-ATG5 ANTIBODY (ABP1038)

Rabbit Polyclonal Antibody

ABP1038

Product Overview

Name	ANTI-ATG5 ANTIBODY (ABP1038)
Catalog #	ABP1038
Clonality	Polyclonal
Accession(Primary)	Q9H1Y0
Other Accession	Q99J83
Application Note (Approx.)	WB1:1000
Precautions	

Target information(Q9H1Y0)

Synonyms


Gene ID 9474

Other Names

Function

Cellular location

Note



Samples: Human (HeLa, MDA-MB-231, K562) cell lysates. PVDF membrane was probed with Anti-ATG5 Antibody (Catalog# PP1038-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.