

# ANTI-ATG5 ANTIBODY (ABP1038)

*Rabbit Polyclonal Antibody*

ABP1038

## Product Overview

<b>Name</b>	ANTI-ATG5 ANTIBODY (ABP1038)
<b>Catalog #</b>	ABP1038
<b>Clonality</b>	Polyclonal
<b>Accession(Primary)</b>	Q9H1Y0
<b>Other Accession</b>	Q99J83
<b>Application Note (Approx.)</b>	WB1:1000
<b>Precautions</b>	

## Target information(Q9H1Y0)

### Synonyms


**Gene ID** 9474

### Other Names

### Function

### Cellular location

### Note



Samples: Human (HeLa, MDA-MB-231, K562) cell lysates. PVDF membrane was probed with Anti-ATG5 Antibody (Catalog# PP1038-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.