

ANTI-BAK ANTIBODY (ABP1034)

Rabbit Polyclonal Antibody


ABP1034

Product Overview

Name	ANTI-BAK ANTIBODY (ABP1034)
Catalog #	ABP1034
Clonality	Polyclonal
Accession(Primary)	Q16611
Application Note (Approx.)	WB1:2000
Precautions	

Target information(Q16611)

Synonyms**Gene ID** 578**Other Names****Function****Cellular location****Note**



Samples: Human (THP-1, Jurkat, MDA-MB-453, 293, RPMI 8226) cell lysates. PVDF membrane was probed with Anti-BAK Antibody (Catalog# PP1034-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.