

# ANTI-BAP31 ANTIBODY (ABM2021)

*Mouse Monoclonal Antibody*


ABM2021

## Product Overview

<b>Name</b>	ANTI-BAP31 ANTIBODY (ABM2021)
<b>Catalog #</b>	ABM2021
<b>Clonality</b>	Monoclonal
<b>Accession(Primary)</b>	P51572
<b>Application Note (Approx.)</b>	WB1:2000
<b>Precautions</b>	

## Target information(P51572)

**Synonyms****Gene ID** 10134**Other Names****Function****Cellular location****Note**



Samples: Human (A431, NCCIT) cell lysates. PVDF membrane was probed with Anti-BAP31 Antibody (Catalog# PM2021-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.