

# ANTI-BETA ACTIN ANTIBODY (ABM2007)

*Mouse Monoclonal Antibody*

ABM2007

## Product Overview

<b>Name</b>	ANTI-BETA ACTIN ANTIBODY (ABM2007)
<b>Catalog #</b>	ABM2007
<b>Clonality</b>	Monoclonal
<b>Accession(Primary)</b>	P60709
<b>Other Accession</b>	P16858;P04797
<b>Application Note (Approx.)</b>	WB1:2000
<b>Precautions</b>	

## Target information(P60709)

### Synonyms


**Gene ID** 60

### Other Names

### Function

### Cellular location

### Note



Samples: Human (A431, HeLa) and Rat (C6) cell lysates.

PVDF membrane was probed with Anti-Beta Actin Antibody (Catalog# PM2007-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody.

Experiment performed under reducing conditions.