

ANTI-Bi-PHOSPHO-ERK1(T202/Y204) ANTIBODY (ABP1047)

Rabbit Polyclonal Antibody


ABP1047

Product Overview

Name	ANTI-Bi-PHOSPHO-ERK1(T202/Y204) ANTIBODY (ABP1047)
Catalog #	ABP1047
Clonality	Polyclonal
Accession(Primary)	P27361
Application Note (Approx.)	WB1:1000
Precautions	

Target information(P27361)

Synonyms**Gene ID** 5595**Other Names****Function****Cellular location****Note**



Samples: Human (HeLa) and Mouse (NIH/3T3) cell lysates. PVDF membrane was probed with Anti-Bi-Phospho-ERK1(T202/Y204) Antibody (Catalog# PP1047-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.