

# ANTI-CALNEXIN ANTIBODY (ABP1088)

*Rabbit Polyclonal Antibody*

ABP1088

## Product Overview

Name	ANTI-CALNEXIN ANTIBODY (ABP1088)
Catalog #	ABP1088
Clonality	Polyclonal
Accession(Primary)	P27824
Other Accession	P35564;P35565
Application Note (Approx.)	WB1:2000
Precautions	

## Target information(P27824)

### Synonyms


Gene ID	821
---------	-----

### Other Names

### Function

### Cellular location

### Note



Samples: Human (PANC-1) cell lysate, and Mouse (brain, liver) and Rat (liver) tissue lysates. PVDF membrane was probed with Anti-Calnexin Antibody (Catalog# PP1088-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.