

ANTI-CALNEXIN ANTIBODY (ABP1088)

Rabbit Polyclonal Antibody

ABP1088

Product Overview

Name	ANTI-CALNEXIN ANTIBODY (ABP1088)
Catalog #	ABP1088
Clonality	Polyclonal
Accession(Primary)	P27824
Other Accession	P35564;P35565
Application Note (Approx.)	WB1:2000
Precautions	

Target information(P27824)

Synonyms


Gene ID 821

Other Names

Function

Cellular location

Note



Samples: Human (PANC-1) cell lysate, and Mouse (brain, liver) and Rat (liver) tissue lysates. PVDF membrane was probed with Anti-Calnexin Antibody (Catalog# PP1088-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.