

ANTI-CAMK2G ANTIBODY (ABP1054)

Rabbit Polyclonal Antibody

ABP1054

Product Overview

Name	ANTI-CAMK2G ANTIBODY (ABP1054)
Catalog #	ABP1054
Clonality	Polyclonal
Accession(Primary)	Q13555
Other Accession	Q923T9;Q923T9
Application Note (Approx.)	WB1:1000
Precautions	

Target information(Q13555)

Synonyms


Gene ID 818

Other Names

Function

Cellular location

Note



Sample: Human (HACAT) cell lysate. PVDF membrane was probed with Anti-CAMK2G Antibody (Catalog# PP1054-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.