

ANTI-CASP9 ANTIBODY (ABP1009)

Rabbit Polyclonal Antibody


ABP1009

Product Overview

Name	ANTI-CASP9 ANTIBODY (ABP1009)
Catalog #	ABP1009
Clonality	Polyclonal
Accession(Primary)	P55211
Other Accession	Q8C3Q9
Application Note (Approx.)	WB1:2000
Precautions	

Target information(P55211)

Synonyms**Gene ID** 842**Other Names****Function****Cellular location****Note**



Sample: Human (A549) cell lysate. PVDF membrane was probed with Anti-CASP9 Antibody (Catalog# PP1009-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under Sample: Human (A549) cell lysate. PVDF membrane was probed with Anti-CASP9 Antibody (Catalog# PP1009-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Sample: Human (A549) cell lysate. PVDF membrane was probed with Anti-CASP9 Antibody (Catalog# PP1009-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.