

# ANTI-CASP9 ANTIBODY (ABP1009)

*Rabbit Polyclonal Antibody*

ABP1009

## Product Overview

Name	ANTI-CASP9 ANTIBODY (ABP1009)
Catalog #	ABP1009
Clonality	Polyclonal
Accession(Primary)	P55211
Other Accession	Q8C3Q9
Application Note (Approx.)	WB1:2000
Precautions	

## Target information(P55211)

### Synonyms


Gene ID 842

### Other Names

### Function

### Cellular location

### Note



Sample: Human (A549) cell lysate. PVDF membrane was probed with Anti-CASP9 Antibody (Catalog# PP1009-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under Sample: Human (A549) cell lysate. PVDF membrane was probed with Anti-CASP9 Antibody (Catalog# PP1009-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Sample: Human (A549) cell lysate. PVDF membrane was probed with Anti-CASP9 Antibody (Catalog# PP1009-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.