

# ANTI-CASPASE-3 ANTIBODY (ABM2026)

*Mouse Monoclonal Antibody*


ABM2026

## Product Overview

<b>Name</b>	ANTI-CASPASE-3 ANTIBODY (ABM2026)
<b>Catalog #</b>	ABM2026
<b>Clonality</b>	Monoclonal
<b>Accession(Primary)</b>	P42574
<b>Application Note (Approx.)</b>	WB1:1000
<b>Precautions</b>	

## Target information(P42574)

**Synonyms****Gene ID** 836**Other Names****Function****Cellular location****Note**



Sample: Human (CCRF-CEM) cell lysate. PVDF membrane was probed with Anti-Caspase-3 Antibody (Catalog# PM2026-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.