

ANTI-CAVEOLIN-2 ANTIBODY (ABP1061)

Rabbit Polyclonal Antibody

ABP1061

Product Overview

Name	ANTI-CAVEOLIN-2 ANTIBODY (ABP1061)
Catalog #	ABP1061
Clonality	Polyclonal
Accession(Primary)	P51636
Application Note (Approx.)	WB1:1000
Precautions	

Target information(P51636)

Synonyms


Gene ID	858
---------	-----

Other Names

Function

Cellular location

Note



Sample: Human (HUVEC) cell lysate. PVDF membrane was probed with Anti-Caveolin-2 Antibody (Catalog# PP1061-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.