

ANTI-CAVEOLIN-3 ANTIBODY (ABP1053)

Rabbit Polyclonal Antibody


ABP1053

Product Overview

Name	ANTI-CAVEOLIN-3 ANTIBODY (ABP1053)
Catalog #	ABP1053
Clonality	Polyclonal
Accession(Primary)	P56539
Application Note (Approx.)	WB1:500
Precautions	

Target information(P56539)

Synonyms**Gene ID** 859**Other Names****Function****Cellular location****Note**



Sample: Human (uterus) tissue lysate. PVDF membrane was probed with Anti-Calveolin-3 Antibody (Catalog# PP1053-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.