

ANTI-CAVEOLIN ANTIBODY (ABP1023)

Rabbit Polyclonal Antibody


ABP1023

Product Overview

Name	ANTI-CAVEOLIN ANTIBODY (ABP1023)
Catalog #	ABP1023
Clonality	Polyclonal
Accession(Primary)	Q03135
Application Note (Approx.)	WB1:1000
Precautions	

Target information(Q03135)

Synonyms**Gene ID** 857**Other Names****Function****Cellular location****Note**



Sample: Human (A2058) cell lysate. PVDF membrane was probed with Anti-Caveolin Antibody (Catalog# PP1023-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.