

ANTI-CD4 ANTIBODY (ABM2023)

Mouse Monoclonal Antibody


ABM2023

Product Overview

| | |
|-----------------------------------|-----------------------------|
| Name | ANTI-CD4 ANTIBODY (ABM2023) |
| Catalog # | ABM2023 |
| Clonality | Monoclonal |
| Accession(Primary) | P01730 |
| Application Note (Approx.) | WB1:1000 |
| Precautions | |

Target information(P01730)

Synonyms**Gene ID** 920**Other Names****Function****Cellular location****Note**



Sample: Human (CCRF-CEM) cell lysate. PVDF membrane was probed with Anti-CD4 Antibody (Catalog# PM2023-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.