

# ANTI-CD44 ANTIBODY (ABM2001)

*Mouse Monoclonal Antibody*


ABM2001

## Product Overview

<b>Name</b>	ANTI-CD44 ANTIBODY (ABM2001)
<b>Catalog #</b>	ABM2001
<b>Clonality</b>	Monoclonal
<b>Accession(Primary)</b>	P16070
<b>Application Note (Approx.)</b>	WB1:2000
<b>Precautions</b>	

## Target information(P16070)

**Synonyms****Gene ID** 960**Other Names****Function****Cellular location****Note**



Samples: Human (HeLa, PC-3, HT-1080) cell lysates.  
PVDF membrane was probed with Anti-CD44 Antibody (Catalog# PM2001-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.