

ANTI-CD68 ANTIBODY (ABM2024)

Mouse Monoclonal Antibody


ABM2024

Product Overview

Name	ANTI-CD68 ANTIBODY (ABM2024)
Catalog #	ABM2024
Clonality	Monoclonal
Accession(Primary)	P34810
Application Note (Approx.)	WB1:1000
Precautions	

Target information(P34810)

Synonyms**Gene ID** 968**Other Names****Function****Cellular location****Note**



Sample: Human (WiDr) cell lysate. PVDF membrane was probed with Anti-CD68 Antibody (Catalog# PM2024-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.