

# ANTI-CD8A ANTIBODY (ABP1010)

*Rabbit Polyclonal Antibody*

ABP1010

## Product Overview

Name	ANTI-CD8A ANTIBODY (ABP1010)
Catalog #	ABP1010
Clonality	Polyclonal
Accession(Primary)	P01732
Application Note (Approx.)	WB1:1000
Precautions	

## Target information(P01732)

### Synonyms


Gene ID 925

### Other Names

### Function

### Cellular location

### Note



Samples: Human (Jurkat, MOLT-4) cell lysates. PVDF membrane was probed with Anti-CD8A Antibody (Catalog# PP1010-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.