

ANTI-CLATHRIN HEAVY CHAIN ANTIBODY (ABP1024)

Rabbit Polyclonal Antibody

ABP1024

Product Overview

| | |
|----------------------------|--|
| Name | ANTI-CLATHRIN HEAVY CHAIN ANTIBODY (ABP1024) |
| Catalog # | ABP1024 |
| Clonality | Polyclonal |
| Accession(Primary) | Q00610 |
| Application Note (Approx.) | WB1:1000 |
| Precautions | |

Target information(Q00610)

Synonyms


Gene ID 1213

Other Names

Function

Cellular location

Note



Sample: Human (HeLa) cell lysate. PVDF membrane was probed with Anti-Clathrin Heavy Chain Antibody (Catalog# PP1024-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.