

# ANTI-CLATHRIN HEAVY CHAIN ANTIBODY (ABP1024)

*Rabbit Polyclonal Antibody*


ABP1024

## Product Overview

<b>Name</b>	ANTI-CLATHRIN HEAVY CHAIN ANTIBODY (ABP1024)
<b>Catalog #</b>	ABP1024
<b>Clonality</b>	Polyclonal
<b>Accession(Primary)</b>	Q00610
<b>Application Note (Approx.)</b>	WB1:1000
<b>Precautions</b>	

## Target information(Q00610)

**Synonyms****Gene ID** 1213**Other Names****Function****Cellular location****Note**



Sample: Human (HeLa) cell lysate. PVDF membrane was probed with Anti-Clathrin Heavy Chain Antibody (Catalog# PP1024-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.