

ANTI-DESMIN ANTIBODY (ABP1100)

Rabbit Polyclonal Antibody


ABP1100

Product Overview

Name	ANTI-DESMIN ANTIBODY (ABP1100)
Catalog #	ABP1100
Clonality	Polyclonal
Accession(Primary)	P17661
Application Note (Approx.)	WB1:1000
Precautions	

Target information(P17661)

Synonyms**Gene ID** 1674**Other Names****Function****Cellular location****Note**



Sample: Human (NCI-H460) cell lysate. PVDF membrane was probed with Anti-Desmin Antibody (Catalog# PP1100-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.