

# ANTI-ERK2 ANTIBODY (ABM2011)

*Mouse Monoclonal Antibody*

ABM2011

## Product Overview

<b>Name</b>	ANTI-ERK2 ANTIBODY (ABM2011)
<b>Catalog #</b>	ABM2011
<b>Clonality</b>	Monoclonal
<b>Accession(Primary)</b>	P28482
<b>Other Accession</b>	P63085;P63086
<b>Application Note (Approx.)</b>	WB1:1000
<b>Precautions</b>	

## Target information(P28482)

### Synonyms


**Gene ID** 5594

### Other Names

### Function

### Cellular location

### Note



Samples: Human (MCF-7, HepG2), Mouse (NIH/3T3, Neuro-2a) and Rat (C6) cell lysates. PVDF membrane was probed with Anti-ERK2 Antibody (Catalog# PM2011-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.