

ANTI-FBL ANTIBODY (ABP1017)

Rabbit Polyclonal Antibody


ABP1017

Product Overview

Name	ANTI-FBL ANTIBODY (ABP1017)
Catalog #	ABP1017
Clonality	Polyclonal
Accession(Primary)	P22087
Application Note (Approx.)	WB1:1000
Precautions	

Target information(P22087)

Synonyms**Gene ID** 2091**Other Names****Function****Cellular location****Note**



Sample: Human (K562) cell lysate. PVDF membrane was probed with Anti-FBL Antibody (Catalog# PP1017-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.