

# ANTI-GAMMA TUBULIN ANTIBODY (ABP1089)

*Rabbit Polyclonal Antibody*

ABP1089

## Product Overview

Name	ANTI-GAMMA TUBULIN ANTIBODY (ABP1089)
Catalog #	ABP1089
Clonality	Polyclonal
Accession(Primary)	Q8VCK3
Other Accession	Q8VCK3
Application Note (Approx.)	WB1:2000
Precautions	

## Target information(Q8VCK3)

### Synonyms


Gene ID	27175
---------	-------

### Other Names

### Function

### Cellular location

### Note



Samples: Mouse (brain, liver) tissue lysates. PVDF membrane was probed with Anti-Gamma Tubulin Antibody (Catalog# PP1089-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.