

ANTI-GAMMA TUBULIN ANTIBODY (ABP1089)

Rabbit Polyclonal Antibody


ABP1089

Product Overview

Name	ANTI-GAMMA TUBULIN ANTIBODY (ABP1089)
Catalog #	ABP1089
Clonality	Polyclonal
Accession(Primary)	Q8VCK3
Other Accession	Q8VCK3
Application Note (Approx.)	WB1:2000
Precautions	

Target information(Q8VCK3)

Synonyms**Gene ID** 27175**Other Names****Function****Cellular location****Note**



Samples: Mouse (brain, liver) tissue lysates. PVDF membrane was probed with Anti-Gamma Tubulin Antibody (Catalog# PP1089-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.