

ANTI-LAMP1 ANTIBODY(ABM2018)

Mouse Monoclonal Antibody


ABM2018

Product Overview

Name	ANTI-LAMP1 ANTIBODY(ABM2018)
Catalog #	ABM2018
Clonality	Monoclonal
Accession(Primary)	P11279
Other Accession	A0A0D9RXY1
Application Note (Approx.)	WB1:4000
Precautions	

Target information(P11279)

Synonyms**Gene ID** 3916**Other Names****Function****Cellular location****Note**



Samples: Human (A431, HeLa, HT-1080, Jurkat) and African Green Monkey (COS-7) cell lysates. PVDF membrane was probed with Anti-LAMP1 Antibody (Catalog# PM2018-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.