

ANTI-LAMP2 ANTIBODY (ABM2008)

Mouse Monoclonal Antibody


ABM2008

Product Overview

Name	ANTI-LAMP2 ANTIBODY (ABM2008)
Catalog #	ABM2008
Clonality	Monoclonal
Accession(Primary)	P13473
Application Note (Approx.)	WB1:4000
Precautions	

Target information(P13473)

Synonyms**Gene ID** 3920**Other Names****Function****Cellular location****Note**



Samples: Human (HeLa, T47D, Jurkat) cell lysates.
PVDF membrane was probed with Anti-LAMP1 Antibody (Catalog# PM2018-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.