

ANTI-LC3 ANTIBODY (ABP1037)

Rabbit Polyclonal Antibody


ABP1037

Product Overview

Name	ANTI-LC3 ANTIBODY (ABP1037)
Catalog #	ABP1037
Clonality	Polyclonal
Accession(Primary)	Q9GZQ8
Other Accession	Q9CQV6;Q62625
Application Note (Approx.)	WB1:2000
Precautions	

Target information(Q9GZQ8)

Synonyms**Gene ID** 81631**Other Names****Function****Cellular location****Note**



Sample: Mouse (NIH-3T3) cell lysate. PVDF membrane was probed with Anti-LC3 Antibody (Catalog# PP1037-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.