

# ANTI-PHOSPHO-ERBB4/HER4 (Y1162) ANTIBODY (ABP1044)

*Rabbit Polyclonal Antibody*

ABP1044

## Product Overview

Name	ANTI-PHOSPHO-ERBB4/HER4 (Y1162) ANTIBODY (ABP1044)
Catalog #	ABP1044
Clonality	Polyclonal
Accession(Primary)	Q15303
Application Note (Approx.)	WB1:2000
Precautions	

## Target information(Q15303)

### Synonyms


Gene ID 2066

### Other Names

### Function

### Cellular location

### Note



Sample: Human (A431) cell lysate. PVDF membrane was probed with Anti-PHOSPHO-ERBB4/HER4 (Y1162) Antibody (Catalog# PP1044-S1) at 1:500 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.