

ANTI-PHOSPHO-PLB (T17) ANTIBODY (ABP1048)

Rabbit Polyclonal Antibody

ABP1048

Product Overview

Name	ANTI-PHOSPHO-PLB (T17) ANTIBODY (ABP1048)
Catalog #	ABP1048
Clonality	Polyclonal
Accession(Primary)	P26678
Other Accession	P61016
Application Note (Approx.)	WB1:1000
Precautions	

Target information(P26678)

Synonyms


Gene ID 5350

Other Names

Function

Cellular location

Note



Sample: Rat (heart) tissue lysate. PVDF membrane was probed with Anti-Phospho-PLB (T17) Antibody (Catalog# PP1048-S1) at 1:1000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.