

ANTI-PHOSPHO-SMAD4 (T277) ANTIBODY (ABP1045)

Rabbit Polyclonal Antibody

ABP1045

Product Overview

Name	ANTI-PHOSPHO-SMAD4 (T277) ANTIBODY (ABP1045)
Catalog #	ABP1045
Clonality	Polyclonal
Accession(Primary)	Q13485
Other Accession	P97471
Application Note (Approx.)	WB1:500
Precautions	

Target information(Q13485)

Synonyms


Gene ID 4089

Other Names

Function

Cellular location

Note



Sample: Mouse (NIH/3T3) cell lysate. PVDF membrane was probed with Anti-PHOSPHO-SMAD4 (T277) Antibody (Catalog# PP1045-S1) at 1:500 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.