

# ANTI-RAB7A ANTIBODY (ABP1025)

*Rabbit Polyclonal Antibody*

ABP1025

## Product Overview

Name	ANTI-RAB7A ANTIBODY (ABP1025)
Catalog #	ABP1025
Clonality	Polyclonal
Accession(Primary)	P51149
Other Accession	P51150;P09527
Application Note (Approx.)	WB1:2000
Precautions	

## Target information(P51149)

### Synonyms


Gene ID	7879
---------	------

### Other Names

### Function

### Cellular location

### Note



Samples: Human (A431, HeLa) and Mouse (C2C12, NIH/3T3) cell lysates. PVDF membrane was probed with Anti-RAB7A Antibody (Catalog# PP1025-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.