

# ANTI-RAB7A ANTIBODY (ABP1025)

*Rabbit Polyclonal Antibody*


ABP1025

## Product Overview

<b>Name</b>	ANTI-RAB7A ANTIBODY (ABP1025)
<b>Catalog #</b>	ABP1025
<b>Clonality</b>	Polyclonal
<b>Accession(Primary)</b>	P51149
<b>Other Accession</b>	P51150;P09527
<b>Application Note (Approx.)</b>	WB1:2000
<b>Precautions</b>	

## Target information(P51149)

**Synonyms****Gene ID** 7879**Other Names****Function****Cellular location****Note**



Samples: Human (A431, HeLa) and Mouse (C2C12, NIH/3T3) cell lysates. PVDF membrane was probed with Anti-RAB7A Antibody (Catalog# PP1025-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.