

# ANTI-TUBULIN BETA-6 ANTIBODY (ABP1084)

*Rabbit Polyclonal Antibody*

ABP1084

## Product Overview

Name	ANTI-TUBULIN BETA-6 ANTIBODY (ABP1084)
Catalog #	ABP1084
Clonality	Polyclonal
Accession(Primary)	Q9BUF5
Other Accession	Q922F4;Q4QQV0
Application Note (Approx.)	WB1:2000
Precautions	

## Target information(Q9BUF5)

### Synonyms


Gene ID 84617

### Other Names

### Function

### Cellular location

### Note



Samples: Mouse (C2C12, NIH/3T3) and Rat (L6) cell lysates, and Mouse (spleen, thymus) tissue lysates. PVDF membrane was probed with Anti-Tubulin Beta-6 Antibody (Catalog# PP1084-S1) at 1:2000 dilution, followed by HRP conjugated secondary antibody. Experiment performed under reducing conditions.