

# Acetyl-Histone H2B (Lys20) Rabbit mAb


ABZ5257

## Product Overview


<b>Name</b>	Acetyl-Histone H2B (Lys20) Rabbit mAb
<b>Catalog #</b>	ABZ5257
<b>Clonality</b>	Monoclonal
<b>Accession(Primary)</b>	P3377
<b>Application Note (Approx.)</b>	WB1:1000-1:5000 IHCP-1:50-1:100
<b>Precautions</b>	

## Target information(P3377)

**Synonyms****Gene ID** 3018**Other Names****Function****Cellular location****Note**



Western blot analysis of Histone H2B (acetyl K20) in HeLa treated with TSA lysates using Acetyl-Histone H2B (Lys20) antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon using Histone H2B (acetyl K20) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.