

Acetyl-Histone H3 (Lys14) Rabbit mAb


ABZ4315

Product Overview

Name	Acetyl-Histone H3 (Lys14) Rabbit mAb
Catalog #	ABZ4315
Clonality	Monoclonal
Accession(Primary)	P6843
Application Note (Approx.)	WB1:500-1:1000
Precautions	

Target information(P6843)

Synonyms**Gene ID** 8350**Other Names****Function****Cellular location****Note**



Immunocytochemistry analysis of Acetyl-Histone H3 (Lys14) (green) in HeLa using Acetyl-Histone H3 (Lys14) antibody, and DAPI (blue).

Western blot analysis of Acetyl-Histone H3 (Lys14) in PC-12 lysates using Acetyl-Histone H3 (Lys14) antibody. Observed band size: 15kDa.



Western blot analysis of Acetyl-Histone H3 (Lys14) in A549, HL-60 lysates using Acetyl-Histone H3 (Lys14) antibody

