

# Acetyl-Histone H3 (Lys4) Rabbit mAb

ABZ4316

## Product Overview

<b>Name</b>	Acetyl-Histone H3 (Lys4) Rabbit mAb
<b>Catalog #</b>	ABZ4316
<b>Clonality</b>	Monoclonal
<b>Accession(Primary)</b>	P6843
<b>Application Note (Approx.)</b>	WB1:1000-1:5000
<b>Precautions</b>	

## Target information(P6843)

**Synonyms****Gene ID** 8350**Other Names****Function****Cellular location****Note**



Western blot analysis of Acetyl-Histone H3 (Lys4) in 3T3 lysates using Acetyl-Histone H3 (Lys4) antibody.

Western blot analysis of Acetyl-Histone H3 (Lys4) in A549, HL-60 lysates using Acetyl-Histone H3 (Lys4) antibody.



Immunocytochemistry analysis of Acetyl-Histone H3 (Lys4) (green) in HeLa using Acetyl-Histone H3 (Lys4) antibody, and DAPI (blue)

