

# Activin-A Antibody

Activin-A, Polyclonal Rabbit Anti-Human Antibody  
ABP0001

## Product Overview

Name Activin-A Antibody

### Description

Activin-A, Polyclonal Rabbit Anti-Human Antibody

### Synonyms

Inhba, Inhibin beta A, FSH releasing protein.

### Introduction

Activins are homodimers or heterodimers of the different  $\alpha$  subunit isoforms, part of the TGF $\beta$  family. Mature Activin A has two 116 amino acids residues  $\alpha$  subunits ( $\alpha$ A- $\alpha$ A). Activin displays an extensive variety of biological activities, including mesoderm induction, neural cell differentiation, bone remodelling, haematopoiesis, and reproductive physiology. Activin takes part in the production and regulation of hormones such as FSH, LH, GnRH and ACTH. Cells that are identified to express Activin A include fibroblasts, endothelial cells, hepatocytes, vascular smooth muscle cells, macrophages, keratinocytes, osteoclasts, bone marrow monocytes, prostatic epithelium, neurons, chondrocytes, osteoblasts, Leydig cells, Sertoli cells, and ovarian granulosa cells.

### Source

E.Coli.

### Physical Appearance

Lyophilized freeze dried powder.

### Formulation

Human Activin-A was lyophilized from a concentrated 1mg/ml protein solution containing 0.1% TFA.

### Stability

Lyophilized Activin-A although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution Activin-A should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

### Purity

Greater than 95% as observed by SDS-PAGE.

### Amino acid sequence

MGLECDGKVN ICCKKQFFVS FKDIGWNDWI IAPSGYHANY CECECPSHIA GTSGSSLSFH STVINHYRMR  
GHSPFANLKS CCVPTKLRPM SMLYYDDGQN IIKKDIQNMI VEECGCS.

### Biological Activity

Biological activity is assessed by the ability to induce cytotoxicity of MPC-11 cells and was found to be 8.8ng/ml corresponding to a specific activity of  $1.1 \times 10^5$  units/mg.

### Solubility

Human INHBA protein should be reconstituted in distilled pyrogen free water to a concentration of 100ug /ml which can then be further diluted to other aqueous solutions.

### Background

**Title:** Research on Activin A Human Recombinant: Molecular Characteristics, Signaling Pathways, Physiological Functions, and Therapeutic Potential  
**Introduction:** Activin A, a member of the transforming growth factor-beta (TGF- $\beta$ ) superfamily, is a multifunctional cytokine that plays a significant role in various biological processes in the human body. Its involvement in diverse physiological and pathological functions has garnered considerable attention in scientific research. This paper aims to provide an overview of Activin A, encompassing its molecular characteristics, signaling pathways, physiological functions, and therapeutic potential. Activin A is encoded by the INHBA gene and is produced as a precursor protein that undergoes post-translational modifications to generate the mature form. The mature Activin A protein consists of two  $\alpha$ -subunits held together by disulfide bonds. These structural features contribute to its functional properties and interactions with specific receptors. Upon binding to its cell surface receptors, Activin A triggers intracellular signaling cascades, leading to various cellular responses. Canonical SMAD-dependent pathway as well as non-SMAD pathways, such as MAPK/ERK, PI3K/Akt, and JNK signaling, are activated by Activin A. The intricate network of signaling pathways enables Activin A to regulate diverse biological processes, including cell proliferation, differentiation, apoptosis, and tissue homeostasis. Activin A exerts its physiological functions in a tissue-specific manner. It plays a critical role in embryonic development, particularly in organogenesis and patterning. Additionally, Activin A is involved in reproductive biology, where it participates in folliculogenesis, spermatogenesis, and hormonal regulation. It also contributes to neural development, immune system modulation, and skeletal homeostasis. The multifunctional properties of Activin A have positioned it as a potential therapeutic target for various diseases. Its involvement in cancer, neurodegenerative disorders, fibrosis, and reproductive disorders has prompted extensive research to explore its therapeutic potential. Understanding the molecular mechanisms underlying Activin A's actions provides valuable insights for developing innovative therapeutic strategies. In conclusion, Activin A is a versatile cytokine with diverse roles in human biology. This research aims to deepen our understanding of its molecular characteristics, signaling pathways, physiological functions, and therapeutic potential. By elucidating the complexities of Activin A, we strive to pave the way for novel therapeutic interventions in various human diseases.

### Precautions

Activin-A Antibody is for research use only and not for use in diagnostic or therapeutic procedures.