

# Dengue 1 NS1 Antibody

*Polyclonal Rabbit Anti Dengue 1 NS1*  
ABP0006

## Product Overview

Name Dengue 1 NS1 Antibody

### Description

Polyclonal Rabbit Anti Dengue 1 NS1

### Formulation

PBS + 0.02% sodium azide.

### Stability

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

### Purity

Protein is >90% pure.

### Purification Method

Purified by affinity chromatographic technique.

### Precautions

Dengue 1 NS1 Antibody is for research use only and not for use in diagnostic or therapeutic procedures.